

digital
physician scale



model 150-10-5

RICE LAKE[®]
HEALTH SOLUTIONS

'To be the best by every measure'[®]

www.ricelake.com/health

digital physician scale

THE RICE LAKE DIGITAL PHYSICIAN SCALE is designed for both modern appearance and contemporary function, with fresh, clean lines and forward-thinking features for the medical and fitness environments of today and tomorrow.

The Rice Lake Digital Physician Scale prepares you for the future requirements of electronic medical record keeping by letting you easily transmit scale data to external devices. The feature-rich Digital Physician Scale is an incredible value, with premium precision and accuracy in a heavy-duty design.

An integrated height rod tucks neatly away within the scale column, and the large, stable scale base is both comfortable and easily transported by two integrated rubber wheels. Whether powered from an AC adaptor or standard alkaline batteries, this model is easy to use. Its large push-button display lights up for easy viewing and an audible beep sounds at the end of each weigh cycle. Add to that a 550lb capacity and exclusive Motiontrap™ movement compensation technology, and this Rice Lake model is one of the most accommodating physician scales on the market.



150-10-5 Standard Features



- Platform Dimensions: 14.5 in W x 14.5 in L x 3 in H
- Motiontrap™ (movement compensation technology)
- Lb/kg toggle key
- Large LCD display
- 6 “AA” alkaline batteries and AC adaptor included
- Low-battery indication/automatic power-off
- BMI (body mass index) function
- Hold function
- Unit of measure: lb only, kg only, lb/kg
- RS-232 or USB output
- Built-In Height measuring rod measures from: 27.5 in to 82 in x 1/8 in (70 cm to 208 cm x 1 mm)
- Integrated wheels for ease of portability

Warranty

- Two-year limited warranty

Capacity	550lb (250 kg)
Graduations	0.2lb (0.1 kg)

Your Rice Lake Weighing Systems distributor is:

RICE LAKE
HEALTH SOLUTIONS

To be the best by every measure®

230 W. Coleman St. • Rice Lake, WI 54868 • USA
TEL: 715-234-9171 • FAX: 715-234-6967 • www.ricelake.com

An ISO 9001 Registered Company © 2011 Rice Lake Weighing Systems PN 126117 10/11
Specifications subject to change without notice.