

INVACARE®

Insignia™

Lightweight Wheelchair



 **progressive**
medical concepts

5020 b u bowman dr, suite 600 • buford, ga 30518
800.511.2601 • fax 770.831.8193



Yes, you can.®

INVACARE®

Insignia™

Lightweight Wheelchair

The new Invacare® Insignia™ Lightweight Wheelchair will be the lightweight chair of choice for providers and patients alike. Fully adjustable arms, from desk to full length, allow providers to stock one chair instead of two, and standard features include a robust wheel lock, simple adjustable head tube design, sleek caster fork and adjustable angle back. Patients will benefit from having a stylish, lightweight chair that is easy to use and adjusts to their specific needs. The Insignia Wheelchair has it all in one package.



1



Convertible, height adjustable arms

In a single design, the arms offer 4" of height adjustment for maximum flexibility and convert from desk to full length for minimum inventory.

4



Large push button quick release axle

Comfort and convenience.

2



Angle adjustable caster journal

The simplicity of this design lends to no loose parts with a 1-2-3 quick and easy adjustment.

- 1) simply loosen the two inner fasteners,
- 2) adjust the two bottom fasteners using the rear edge of the caster fork as your guide,
- 3) re-tighten the inner fasteners.

3



Robust wheel lock

User, therapist and provider preferred, the robust wheel lock is easy to locate and engage with a "click" that confirms it is locked. It is also easy to maintain.



Standard Features

- Lightweight, aluminum frame
- Quick release axles
- Swingaway footrests
- Push-button, height adjustable arms that offer desk-to-full length changeability
- Height and angle adjustable back includes spring-loaded, flip-up anti-tippers
- Adjustable angle head tubes
- Stock models for quick ship convenience

5



Adjustable angle back

Industry proven bolt-through bracket adjusts from -5° to 20° in 5° increments.

6



Flip-back arm features a large paddle release

Ergonomically correct, the large paddle release is easy to operate and requires minimum dexterity and is removable.

7



Spring-loaded, flip-up anti-tippers

Allow the caregiver to pull and flip without removing for easy navigation over curbs and thresholds.

8



Height adjustable back

The back features a one piece, tool-free height adjustment knob that offers adjustment from 15" – 19" in 1" increments.



Insignia™ Wheelchair Specifications

Overall	
Width Open	7.50" plus seat width
Width Closed	12"
Depth (w/o Riggings)	31"
Seat Widths	16", 18", 20"
Seat Depths	16", 18"

Knee-to-Heel* (range adjusts in 1" increments)

	Min.	Max.
60°	16"	20"
Elevating	17"	20.5"

* Measurement taken from top of seat rail to back of footplate

Arm Type**

Space-Saver Dual-Point Desk-to-Full Length, Adjustable Height	Min.	Max.
	9"	13"

**Measurement taken from top of seat rail to top of arm pad

Back Heights	15" – 19" in 1" increments
Weight Capacity	250 lb.
Seat-to-Floor	19.5", 18.5", 17.5" or 16.5" with 22"/6" wheel kit

DISTRIBUTED BY



5020 b u bowman dr, suite 600 • buford, ga 30518
800.511.2601 • fax 770.831.8193

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

© 2010

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted.
Specifications are subject
to change without
notification.

Form No. 09-101
Rev. 4/10

