

A cost-effective chair, the LT offers streamlined features focusing on BRODA's core elements:


- BRODA's Comfort Tension Seating® aids in long term comfort and pressure ulcer prevention.
- Unlike traditional recliners, the LT offers smooth infinitely adjustable seat tilt for individualized positioning to help prevent falls and redistribute pressure.
- Unique front pivot seat tilt enables patients to look forward, even in full tilt, for socializing with others - Enhancing quality of life.
- Lightweight, yet exceptionally durable design provides ease of maneuverability.



28° of infinitely adjustable seat tilt provides safe and comfortable positioning.



Versatility allows patients to sit up at a table or rest in the chair.



*Providing peace of mind through unsurpassed comfort and pressure redistribution*

## The Ideal Chair for Hospice Care and General Patient Seating

### STANDARD FEATURES

- BRODA Comfort Tension Seating®
- 16 gauge powder coated tubular steel
- 19" seat width
- 18" seat depth
- 32" back height
- 28° infinitely adjustable seat tilt
- 20" seat to floor height
- Infinitely adjustable elevating legrest
- Soft, durable polyurethane armrests
- 3 x 5" total lock casters, 1 x 5" rear directional Lock caster
- 26" overall width
- 65 lbs. overall weight
- 250 lbs. weight capacity

### STANDARD REMOVABLE VINYL PADDING:

- Neckrest
- Side panels
- Seat and back pads
- Legrest pad



Manufactured in North America where we produce all our products, the LT offers true BRODA quality and safety.



Legrest is infinitely adjustable for individual positioning and additional pressure redistribution.



*The LT helps you to provide the best quality of care for patients:*

Removable vinyl padding is easy to clean.

Caster configuration allows for smooth tracking and ease of steering in tight corners.

Small footprint compared to traditional recliners, is great in small spaces.

Durable frame and components allow for overall ease of use.

*The LT offers exceptional positioning options that provide comfort, stability and safety to help reduce slumping, sliding and falls.*

### BRODA Comfort Tension Seating®

The Key to Comfort is in BRODA's Strapping! Our strapping contours to the body helping to prevent skin breakdown and provide long term sitting comfort.

- Creates air flow for heat and moisture exchange
- Memory retention allows straps to return to their original shape within seconds.

