

SCALES FOR HEALTH AND FITNESS



Service. Selection. Speed.

Rice Lake Weighing Systems



Rice Lake Weighing Systems, an [ISO 9001](#)-registered company, and a leader in delivering high quality weighing and process-control solutions to industry for 60 years. We design, manufacture, and distribute a wide variety of weighing scales and processing control systems, and offer the largest selection of load cells, weigh modules, and advanced instrumentation---over 20,000 products---to meet virtually every weighing requirement.

We also manufacture precision calibration weights, including all ASTM, OIML, and NBS classes for balances and scales. Traceable and Calibration Certificates are available with every weight we sell. Our state-of-the-art metrology lab is accredited by the National Voluntary Laboratory Accreditation Program (NVLAP) and provides quality weight services, including recalibration and refurbishing of existing weights.

Our newly designed Medical and Fitness Scale line is ready for launch, and offers features and functionality that healthcare and fitness professionals have come to expect.

Our commitment to our customers is how **Rice Lake Weighing Systems** has built its reputation, and we are confident we can bring our core capabilities to healthcare and fitness..beyond measurement!

Digital Physician Scales

Digital Stand On Scales offer faster, more accurate weighing. High quality, feature rich scales are designed to work in all medical and fitness applications. Increase productivity and reduce costs and down-time.

Eye Level Digital Stand On Scale

Model Number 140-10-5 Part Number 100124

- 550 lb/250 kg Capacity x .2 lb/.1 kg Resolution
- Large 14.6" x 15" x 3" Base Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- BMI-Weight Hold-lb/kg Switch/Keypad Tare
- RS-232 Output to EMR or Software
- Height Rod 35" to 82" x 1'8" with metric
- Program Weight Lock Out Mode
- Updateable Firmware
- Wheels for Mobility
- 2 Year Warranty



Digital Physician Scales

Digital Waist High Scale offers faster, more accurate weighing. High quality, feature rich scales are designed to work in all medical and fitness applications. Increase productivity and reduce costs and down-time.

Waist-High Digital Stand On Scale

Model Number 140-10-6 Part Number 100189

- 550 lb/250 kg Capacity x .2 lb/.1 kg Resolution
- Large 14.6" x 15" x 3" Base Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- BMI-Weight Hold-lb/kg Switch/Keypad Tare
- RS-232 Output to EMR or Software
- Program Weight Lock Out Mode
- Updateable Firmware
- Wheels for Mobility
- 2 Year Warranty



Digital Physician Scales

Remote Display Scale offers faster, more accurate weighing. High quality, feature rich scales are designed to work in all medical and fitness applications. Increase productivity and reduce costs and down-time.

Remote Display Digital Stand On Scale

Model Number 140-10-7 Part Number 100185

- 550 lb/250 kg Capacity x .2 lb/.1 kg Resolution
- Large 14.6" x 15" x 3" Base Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- BMI-Weight Hold-lb/kg Switch/Keypad Tare
- RS-232 Output to EMR or Software
- Program Weight Lock Out Mode
- Updateable Firmware
- 2 Year Warranty



Balance Beam Scale

Mechanical Balance Beam Scale for accurate weighing. High quality with steel base and cast iron weigh bridge for long lasting results. Dual reading in pounds and kilograms includes dual reading Height Rod and Wheels for mobility.

Balance Beam

Model Number RL-MPS Part Number 102613

- 440 lb/200 kg Capacity x 1/4 lb/100 g Resolution
- Large 10.5" x 14" x 3" Base Size
- Height Rod 24" to 84" x 1/8" and metric equivalent
- Wheels for Mobility
- 2 Year Warranty



Digital Pediatric Scales

Digital accuracy offers faster, more reliable weighing. High quality, feature rich scales are designed to work in all medical applications. Increase productivity and reduce costs and down-time.

Digital Pediatric Scale

Model Number RL-DBS Part Number 101133

- 44 lb/20 kg Capacity x .5 oz/.10 g Resolution
- Large 23.5" x 11" Tray Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- Weight Lock-lb/kg Switch/Keypad Tare



Handrail Scale

Digital Handrail Scale offers faster, more accurate weighing. High quality, feature rich scales are designed to work in all medical and fitness applications. Increase productivity and reduce costs and down-time.

Handrail Scale

Model Number 240-10-7 Part Number 100193

- 700 lb/310 kg Capacity x .2 lb/.1 kg
- Large 20" x 25" x 3" Base Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- RS-232 Output
- Program Weight Lock Out Mode
- Updateable Firmware
- Collapsible Rails for Storage
- Wheels for Mobility
- 2 Year Warranty



Folding Wheelchair Scale

Our new Folding Wheelchair Scale offers faster, more accurate weighing. High quality, feature rich scales with large base and 800 lb. capacity. Designed to work in all medical and fitness applications. Increase productivity and reduce costs and down-time.

Folding Wheelchair Scale

Model Number 340-10-2 Part Number 111712

- 1000 lb/500 kg Capacity x .2 lb/.1 kg Resolution
- Large 29.5" x 32" x 2.5" Base Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- RS-232 Output
- Keypad Tare for Wheelchair
- Program Weight Lock Out Mode
- Updateable Firmware
- Folding Platform
- Wheels for Mobility
- 2 Year Warranty



Wheelchair Ramp Scale

New Wheelchair Ramp Scale offers faster, more accurate weighing. High quality, feature rich scales with large base and **800 lb.** capacity. Designed to work in all medical and fitness applications. Increase productivity and reduce costs and down-time.

Wheelchair Ramp Scale

Model Number 340-10-7 PN 111713

- 800 lb/400 kg Capacity x .2 lb/.1 kg Resolution
- Large 29.5" x 32" x 2.5" Base Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- RS-232 Output
- Keypad Tare for Wheelchair
- Program Weight Lock Out Mode
- Updateable Firmware
- Wheels for Mobility
- 2 Year Warranty



Bariatric Wheelchair Ramp Scale

Bariatric Wheelchair Ramp Scale offers faster, more accurate weighing with dual ramps. High quality, feature rich scales with large base and 1000 lb. capacity. Designed to work in all medical and fitness applications. Increase productivity and reduce costs and down-time.

Bariatric Wheelchair Ramp Scale

Model Number SB-1000 Part Number 106658

- 1000 lb/500 kg Capacity x .2 lb/.1 kg Resolution
- Large 36" x 36" x 2.5" Platform Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- RS-232 Output
- Keypad Tare for Wheelchair
- Program Weight Lock Out Mode
- Updateable Firmware
- 2 Year Warranty



Dialysis Wheelchair Ramp Scale

Dialysis Wheelchair Ramp Scale offers faster, more accurate weighing built in to the floor. High quality, feature rich scales with large base and 1000 lb. capacity. Designed to work in all medical and fitness applications. Increase productivity and reduce costs and down-time.

Dialysis Wheelchair Ramp Scale

Model Number SD-1000 Part Number 106789

- 1000 lb/500 kg Capacity x .2 lb/.1 kg Resolution
- Large 36" x 36" x 2.5" Platform Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- RS-232 Output
- Keypad Tare for Wheelchair
- Program Weight Lock Out Mode
- Updateable Firmware
- 2 Year Warranty



Chair Scales

All New Digital Chair Scales for faster, more accurate weighing. High quality, feature rich scales that help increase productivity and reduce costs and down-time. High Capacity **600 lb.** scale is accurate and dependable.

Digital Chair Scale

Model Number 540-10-1 Part Number 100232

- 600 lb/270 kg Capacity x .2 lb/.1 kg Resolution
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- RS-232 Output
- Program Weight Lock Out Mode
- Updateable Firmware
- Wheels for Mobility
- 2 Year Warranty

Mechanical Chair Scale

Model RL-MCS Part Number 111326

- 440 lb/250 kg Capacity x .2 lb/.1 kg Resolution
- Beam reading front and back side
- Movable arm rests
- Wheels for Mobility
- 2 Year Warranty



Premium Chair Scale

All New Premium Digital Chair Scales for faster, more accurate weighing. High quality, feature rich scales that help increase productivity and reduce costs and down-time. High Capacity **600 lb.** scale is accurate and dependable.

Digital Chair Scale

Model Number 540-10-2 Part Number 101664

- 600 lb/270 kg Capacity x .2 lb/.1 kg Resolution
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- RS-232 Output
- Program Weight Lock Out Mode
- Updateable Firmware
- Wheels for Mobility
- Padded Seat with waterproof cover
- 2 Year Warranty



Neonatal Scale

Fine resolution Neonatal Scale for faster, more accurate weighing. Accurate and durable and designed to meet the strictest standards in healthcare. Increase productivity and reduce costs and down-time.

Neonatal Scale

Model Number 610-10-2 Part Number 101133

- 30 lb/15 kg Capacity x .1 oz/.1g Resolution
- Large 12" x 22" x 4 1/4" Tray Size
- A/C and Battery Power-Adapter Included
- Measuring Tape on Tray
- Motion Compensation Technology
- RS-232 Output
- 2 Year Warranty



Platform/Veterinary Scale

Dual use Platform/Vet Scale for faster, more accurate weighing. High Quality, Feature Rich Scales that help increase productivity and reduce costs and down-time.

Platform/Vet Scale

Model Number 740-10-2 Part Number 100643

- 600 lb/270 kg Capacity x .2 lb/.1 kg Resolution
- Large 22" x 42" x 2" Platform Size
- 1" LCD Display
- A/C and Battery Power-Adapter Included
- RS-232 Output
- Program Weight Lock Out Mode
- Updateable Firmware
- Wheels for Mobility
- 2 Year Warranty

