



 **progressive**
medical concepts



The **DYN**amic solution for safer mobility,
ERGOnomic positioning, and improved comfort.

Provides better positioning and pressure management than ordinary wheelchairs.

Reduces the risk of falls and facilitates independent mobility.

Fully adjusts for a perfect fit.

**SERIES -S-
MOBILITY
CHAIR**

DYN-ERGO Scoot Chair

**SERIES -S-
MOBILITY
CHAIR**

DYN-ERGO Scoot Chair

The **DYN**amic solution for safer mobility,
ERGOnommic positioning, and improved comfort.

SAFETY, MOBILITY, POSITIONING, AND COMFORT.

Designed to reduce the risk of falls for wheelchair users who propel themselves with their feet, the **DYN-ERGO Scoot Chair** safely and comfortably facilitates independent mobility. A complete range of adjustments including armrest height, seat depth, recline angle, and tilt angle make it easy to achieve proper positioning while providing comfortable and safe independent mobility.

STANDARD FEATURE :

- Adjustable seat height.
- Adjustable seat depth.
- Adjustable recline angle seat to back.
- Adjustable armrest height.
- Adjustable tilt angle.

LARGE SIZE AVAILABLE



AVAILABLE OPTIONS & ACCESSORIES :

- SA-10** - Self-storing and Removable Footrests.
- SA-20** - Structural Lateral Supports adjust to provide comfortably correct trunk positioning.
- SA-30** - 13" rear wheels with handbrakes.
- SA-35** - Steel Caster option. Includes 2 rear total lock swivel casters, two front swivel casters - one with directional lock.
- SA-40** - Removable Activity Tray.
- SA-50** - Adjustable Cervical Headrest.
- SA-60** - PRE-CURE cushion allows the insertion of specially pressure management cushions. (See PRE-CURE's brochure)
- SA-70** - Heel Support Loops.
- SA-80** - Increased seat depth option 15.5" to 20.75" (39.4 cm to 52.7 cm).
- SA-90** - Increased seat height option 16.5" to 21.5" (41.9 cm to 54.6 cm).



TECHNICAL SPECIFICATIONS S - 150

- Max. length - reclined and tilted w/footrest extended :
- Min. length - back upright w/o footrest :
- Max. width w/13" wheels and footrest option :
- Min. width w/o accessories :
- Max. height w/back upright and seat height at 18.5" :
- Min. height w/back reclined and seat height at 13.5" :
- Seat width between armrests :
- Seat depth :
- Seat height :
- Armrest height from seat :
- Seat-to-backrest angle :
- Max. tilt (depends on adjustments):
- Weight of chair w/o accessories :
- Weight capacity :
- Ships :



Standard Width



24" Width

58" / 147.3 cm	58" / 147.3 cm
31" / 78.7 cm	31" / 78.7 cm
28" / 71.1 cm	33" / 83.8 cm
26" / 66 cm	31" / 78.7 cm
51.5" / 130.8 cm	51.5" / 130.8 cm
37" / 94 cm	37" / 94 cm
18.5" / 47 cm	23.5" / 59.7 cm
Adjusts from 13.75" to 18" / 34.9 cm to 45.7 cm	
Adjusts from 13.5" to 18.5" / 34.3 cm to 47 cm	
Adjusts from 8" to 9.5" / 20.3 cm to 24.1 cm	
95° to 115°	95° to 115°
34°	34°
49 lbs / 22.3 kg	65 lbs / 29.5 kg
250 lbs / 114 kg	400 lbs / 182 kg
KD via UPS	Common Carrier

Contract grade leatherette upholstery standard. Crypton™ fabric available.
Conforms to CA Technical Bulletin 117.

5020 B U Bowman Dr., Suite 600
Buford GA 30518
sales@progressivemedical
concepts.com

800-511-2601

