

Florien II

A dynamic and highly mobile seating system, the Florien II enables greater independence by promoting correct posture and facilitating a supportive, comfortable seated position.

Unique features including an ergonomically moulded backrest, removable arms and modern styling make this an excellent solution for a wide variety of users and environments including long term care facilities.

General Use

- Long term care facilities
- General care of the elderly wards and bedside use
- Domestic home environments
- Hospices
- For those that side transfer
- General nursing care for the semi and non-ambulant user



Features & Benefits

- Manual Tilt-in-Space positioning
- Intelli-Gel® Integral pressure relieving seat module
- Removable Arms for easier side transfers
- Adjustable Lateral Supports
- Manual Legrest elevation
- Flip Up Footrest
- 25 degree negative Legrest angle to facilitate simpler standing transfers
- Three seat width options: 16", 18" or 20" (determined by width of arm module)
- Articulated integral head support
- Supremely comfortable moulded backrest
- Multi-directional stretch, breathable Dartex material on all user contact surfaces
- Ergonomically positioned push handle
- One directional rear lock castor and three lockable castors.

Weight Limit

- 390 lbs

Overall dimensions:

- 35" with legrest (l) x 28" (w) x 52" (h)

Approximate product weight:

- 121 lbs