




**This Entry-Level LTC Mattress has a Visco Foam Heel Section and a 7-Year Warranty!**



**COMPARE TO**  
 \*Direct Supply® Panacea® Tender  
 Invacare® Solace Prevention  
 Span America® GeoMatt Plus Blue  
 Chip® Relief Care Pro

## Gravity 7 Long Term Care Pressure Redistribution Mattress

Item #	Raised Side Rails	Product Description	Product Dimensions	Product HCPCS	UOM
15770	No	Multi-Layered/Multi-zoned Foam Mattress	80" (L) x 36" (W) x 6" (H)	E0184	1/ea
15777	Yes	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	80" (L) x 36" (W) x 6" (H)	E0184	1/ea
15784	No	Multi-Layered/Multi-zoned Foam Mattress	84" (L) x 36" (W) x 6" (H)	E0184	1/ea
15785	Yes	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	84" (L) x 36" (W) x 6" (H)	E0184	1/ea
15876	No	Multi-Layered/Multi-zoned Foam Mattress	76" (L) x 36" (W) x 6" (H)	E0184	1/ea
15877	Yes	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out	76" (L) x 36" (W) x 6" (H)	E0184	1/ea

\*Direct Supply® and Panacea® are registered trademarks of Direct Supply, Inc.

specifications	
Cover Material	Fluid-Resistant Nylon Top   Waterproof Vinyl Bottom
Mattress Layers of Foam	3
Side Rails	Yes
Dedicated Heel Section	Yes
Fire Rating	CFR 16 Part 1633 Federal Fire Code
Weight Capacity	300 lbs.
Warranty	7 Year Limited Foam   1 Year Cover

- Economical option maximizes value and comfort with 3 full layers of foam available with or without raised side rails.
- Top layer uses unique die-cutting and a true memory foam heel section to optimize zoned pressure redistribution.
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress significantly.
- Cover is fire retardant, fluid proof, tear resistant, low shear, and easy to clean.