

NEW & IMPROVED!

VIPER PLUS GT

ALL THE FEATURES
OF THE ORIGINAL
VIPER PLUS GT

PLUS

THE NEW "BUILT IN"
SEAT EXTENSION!

**MOST
VERSATILE
CHAIR OUT OF
THE CARTON!**



Patent Pending
Aluminum Viper Plus GT
Deluxe High Strength, Lightweight,
Dual Axle (K0004)

FEATURES:

- Adjustable seat rail extension and upholstery with 2" extension easily adjusts from 16" to 18" depth.
- **The Viper Plus GT offers the additional convenience of having the ("Rapid Delivery Program-RDP"), all in one carton and ready to go.**
- All Aluminum frame.
- Dual Axle allows for multiple seat to floor height positions.
- Adjustable height, flip back, removable, padded arms.
- Weighs under 32 lbs. (Excluding front rigging).
- Over center cross brace eliminates seat guides, works with seating and positioning accessories.
- Adjustable angle back easily adjusts from 5° - 20°.
- Adjustable height back adjusts from 17" - 19" in 1" increments, easily works adaptable for individual patient needs.
- Durable, heavy gauge flame retardant nylon upholstery with velcro modesty flaps, resists mildew and bacteria.
- Quick release, 24" composite Mag-style wheels are lightweight and maintenance free.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Adjustable angle caster forks have 3 height adjustments and angle adjustability.
- Tool free adjustable swing away footrests also feature easy to access release handles.
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight.
- Adjustable anti-tipper with flip-up wheels. Never needs to be removed.
- Comes with push to lock wheel locks.
- Patent Pending.
- Limited Lifetime Warranty.

NEW!

Step 1



Step 2



Step 3



**SEAT RAIL TELESCOPES
AND UPHOLSTERY
EXTENDS EASILY
FROM 16" TO 18"!**

wheelchairs
To order call toll free: 877.224.0946

drive[™]

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*

Viper Plus GT – DELUXE HIGH STRENGTH, LIGHTWEIGHT, DUAL AXLE

| | | | |
|--------------------|--|--------------------------------------|------|
| PLA416FBFAARAD-SF | 16" Flip Back, Detachable, Adj. Height Full Arm, Swing-away Footrest | (comes with built-in seat extension) | 1/cs |
| PLA416FBFAARAD-ELR | 16" Flip Back, Detachable, Adj. Height Full Arm, Elevating Legrest | (comes with built-in seat extension) | 1/cs |
| PLA416FBDAARAD-SF | 16" Flip Back, Detachable, Adj. Height Desk Arm, Swing-away Footrest | (comes with built-in seat extension) | 1/cs |
| PLA416FBDAARAD-ELR | 16" Flip Back, Detachable, Adj. Height Desk Arm, Elevating Legrest | (comes with built-in seat extension) | 1/cs |
| PLA418FBDAARAD-SF | 18" Flip Back, Detachable, Adj. Height Desk Arm, Swing-away Footrest | (comes with built-in seat extension) | 1/cs |
| PLA418FBDAARAD-ELR | 18" Flip Back, Detachable, Adj. Height Desk Arm, Elevating Legrest | (comes with built-in seat extension) | 1/cs |
| PLA418FBFAARAD-SF | 18" Flip Back, Detachable, Adj. Height Full Arm, Swing-away Footrest | (comes with built-in seat extension) | 1/cs |
| PLA418FBFAARAD-ELR | 18" Flip Back, Detachable, Adj. Height Full Arm, Elevating Legrest | (comes with built-in seat extension) | 1/cs |
| PLA420FBDAARAD-SF | 20" Flip Back, Detachable, Adj. Height Desk Arm, Swing-away Footrest | (comes with built-in seat extension) | 1/cs |
| PLA420FBDAARAD-ELR | 20" Flip Back, Detachable, Adj. Height Desk Arm, Elevating Legrest | (comes with built-in seat extension) | 1/cs |
| PLA420FBFAARAD-SF | 20" Flip Back, Detachable, Adj. Height Full Arm, Swing-away Footrest | (comes with built-in seat extension) | 1/cs |
| PLA420FBFAARAD-ELR | 20" Flip Back, Detachable, Adj. Height Full Arm, Elevating Legrest | (comes with built-in seat extension) | 1/cs |
| PLA422FBDAAR-SF | 22" Flip Back, Detachable, Adj. Height Desk Arm, Swing-away Footrest | | 1/cs |
| PLA422FBDAAR-ELR | 22" Flip Back, Detachable, Adj. Height Desk Arm, Elevating Legrest | | 1/cs |
| PLA422FBFAAR-SF | 22" Flip Back, Detachable, Adj. Height Full Arm, Swing-away Footrest | | 1/cs |
| PLA422FBFAAR-ELR | 22" Flip Back, Detachable, Adj. Height Full Arm, Elevating Legrest | | 1/cs |

REPLACEMENT FRONT RIGGINGS (All Viper Wheelchairs Include Front Riggings.)

| | | |
|---------|---|--------|
| AL3-4SF | Swing-away Footrest (17"- 21") Footplate Extensions with Heel-loops | 1pr/cs |
| ALELR | Swing-away Elevating Legrest (17.5"- 21") with Padded Calf Pads | 1pr/cs |

ARMS (Traditional Arms style for Viper Plus GT)

| | | |
|------------|--|------|
| STDS1A35DR | Right, Detachable Flip Back, Desk Arm, Adjustable Height | 1/ea |
| STDS1A35DL | Left, Detachable Flip Back, Desk Arm, Adjustable Height | 1/ea |
| STDS1A35FR | Right, Detachable Flip Back, Full Arm, Adjustable Height | 1/ea |
| STDS1A35FL | Left, Detachable Flip Back, Full Arm, Adjustable Height | 1/ea |

WHEEL LOCKS (Push to Lock)

| | | |
|------------|---|------|
| STDS1A351R | Right, Push to Lock, Detachable and Flip Back Arm | 1/ea |
| STDS1A351L | Left, Push to Lock, Detachable and Flip Back Arm | 1/ea |

UNIVERSAL FORWARD STABILIZER

STDS860

ANTI TIPPERS WITH WHEELS

STDS807

SEAT BELTS

STDS850 (Auto Clasp Type)

STDS851 (Velcro Type Closure)

BIARIATRIC SEAT BELTS

STDS855 (Auto Clasp Type)

STDS856 (Velcro Type Closure)

TELESCOPING I.V. POLE ATTACHMENT

STDS820

UNIVERSAL OXYGEN "E" CYLINDER HOLDER WITH I.V. POLE

STDS804

6" WHEEL LOCK EXTENSIONS

STDS801GT

OBLIQUE PROJECTION HANDRIM

STDS1022

LIGHTWEIGHT 8" CASTER

STDS1087

FLAT FREE WHEEL WITH HANDRIM

STDS840

LIGHTWEIGHT 24" WHEEL

STDS1085

STDS1617N (For use with 17" Viper Plus GT)

STDS1618N (For use with 18" Viper Plus GT)

STDS1619N (For use with 19" Viper Plus GT)

OVERHEAD ANTI THEFT DEVICE

STDS823

STDS821 (With I.V. Hooks)

STDS834 (Single Pole)

WHEELCHAIR CANE/CRUTCH HOLDER

STDS1034

HEEL LOOPS

STDS831

LIMB SUPPORT

WASR Right

WASL Left

CHART CARRY POCKET

STDS835 (For use with all 16", 18" and 20" Wheelchairs)

SUPER HEMI KIT

STDS5HK

ONE ARM DRIVE

STDS4055R, STDS4055L (For use with 18" Wheelchairs Only)

SPECIFICATIONS

**Adjustable on new model #'s. Available in 2011.

| Seat Width | Overall Width (open) | Overall Width (closed) | Overall Height | Seat To-Floor Height | Back Height | Width Seat | Depth Seat **(adjustable) | Overall Length w/Riggings | Product Weight | Carton Shipping Weight | Weight Limit |
|------------|----------------------|------------------------|----------------|----------------------|-------------|------------|------------------------------|---------------------------|----------------|------------------------|--------------|
| 16" | 24" | 12.5" | 36" | 17.5"-19.5" | 17"-19" | 16" | 16"-18" | 42" | 31.4 lbs. | 36.4 lbs. | 300 lbs. |
| 18" | 26" | 12.5" | 36" | 17.5"-19.5" | 17"-19" | 18" | 16"-18" | 42" | 31.9 lbs. | 36.9 lbs. | 300 lbs. |
| 20" | 28" | 12.5" | 36" | 17.5"-19.5" | 17"-19" | 20" | 16"-18" | 42" | 33.4 lbs. | 38.4 lbs. | 300 lbs. |
| 22" | 30" | 12.5" | 36" | 17.5"-19.5" | 17"-19" | 22" | 18" | 42" | 39.5 lbs. | 44.5 lbs. | 300 lbs. |

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

wheelchairs
To order call toll free: 877.224.0946

drive™